

STATCLEAN TECHNOLOGY (S) PTE. LTD.

30 Kaki Bukit Road 3, # 03-12/13/14, Empire Technocentre, Singapore 417819. Tel: 65-6746 2706 Fax: 65-6746 3045 / 65-6746 3755 Email: sales@statclean.com Website: www.statclean.com

CPE Shoe Cover (STDP-SC/CPE/1540)

Specifications:

Material : 91% Low Denisty Polythene Film(Embossed) + 9% Colour

Masterbatch - Food Contact Quality

Enclosed Latex Elastic Thickness – 0.025mm

Colour : Blue

Storage : Temperature - Between -5°C and 40°C

Humidity - Below 80%

Sunlight - Not to expose to direct sunlight

Frequency - To be kept away from high frequency equipment

Shelf Life : 3-5 years if store as above

Sterilization : Gamma Radiation - Possible – not advised

Causes colour change and increased ageing

Autoclave - No ETO - Possible

> Gas Contact – 3 hours/850mg Evaporation 96 hours at 40°C

Intended Use : General Medical and Non Medical Use.

Food Environments

Classification : CE ACC. MPF Class 1 Product

DIN EN ISO 9002 By Importer

Risk Analysis : Inflammable, Strictly avoid contact with open fire and heat, only handle

non toxic and not aggressive/dangerous materials

Packaging: 100pcs per poly bag, 20 poly bags per case

Dimension : A - Length : 400mm (+/-5%)

B - Width : 150mm (+/- 5%)

A