

## STATCLEAN TECHNOLOGY (S) PTE. LTD.

30 Kaki Bukit Road 3, # 03-12/13/14, Empire Technocentre, Singapore 417819. Tel: 65-6746 2706 Fax: 65-6746 3045 / 65-6746 3755 Email: sales@statclean.com Website: www.statclean.com

## LDPE Oversleeves (Elasticated at cuff and elbow) (STDP-SLEEVEC/2040)

**Specifications:** 

Material : 91% Low Denisty Polythene Film + 9 % Colour

Masterbatch – Food Contact Quality

Enclosed Latex Elastic Thickness – 0.02mm

Colour : Blue

Storage : Temperature - Between -5°C and 40°C

Humidity - Below 80%

Sunlight - Not to expose to direct sunlight

Frequency - To be kept away from high frequency equipment

**Shelf Life** : 3-5 years if store as above

**Sterilization** : Gamma Radiation - Possible – not advised

Causes colour change and increased ageing

Autoclave - No ETO - Possible

> Gas Contact – 3 hours/850mg Evaporation 96 hours at 40°C

**Intended Use** : General Medical and Non Medical Use.

Food Environments

Classification : CE ACC. MPF Class 1 Product

DIN EN ISO 9002 By Importer

**Risk Analysis** : Inflammable, Strictly avoid contact with open fire and heat, only handle

non toxic and not aggressive/dangerous materials

**Packaging** : 100 pieces per poly bag, 20 poly bags per case

OR

50 pieces per poly bag (w/o elastic bands), 40 poly bags per case

 $\begin{tabular}{ll} \textbf{Dimension} & : A-Length & : 400mm (+/-5\%) \\ \end{tabular}$ 

B – Width : 200 mm (+/-5%)

A B